Solid State Refractive Laser OLIMP-2000/213

РЫБИНСКИЙ ЗАВОД ПРИБОРОСТРОЕНИЯ







Specifications

Source of radiation: solid state pulse Nd: YAG with diode pumping

	and Q-SW modulation	
Wavelength:	213 nm	
Pulse duration:	10 nsec	
Pulse generation frequency:	300 Hz	
Transverse profile of energy distribution:	Gaussian, TEMoo	
Laser beam in focus diameter:	0.5 mm	
The maximum output laser energy:	1 mJ	
Forming system:	"flying" spot	
Auto-tracking:	on the limb and iris	
Ablation zone diameter:	to 20 mm	
Reaching operating mode:	30 min	
Time of continuous operation:	12 hours	
Weight:	160 kg	
External dimensions (l×w×h):	800 mm×1000 mm×1200 mm	
Laser cooling system:	liquid-filled	
Power supply voltage:	220 V, 50 Hz	
Power consumption:	0.9 kWh	
Service life:	at least 10 years	

Clinical possibilities

Myopia up to	-15 diopters	
Hypermetropia up to	+6 diopters	
Myopic astigmatism up to	-6 diopters	
Hypermetropic astigmatism up to	+6 diopters	
Glaucoma surgery		

Types of surgery made:

LASIK, PRK, MAGEK, Epi-LASIK, LASEK, PTK Personalized ablation

Registration Certificate No. ΦCP 2010/08230 Certificate of conformity № POCC RU.ИM25.B03771







Solid state technology

A solid state nanosecond pulsed Nd:YAG laser is used for performing refractive operations with the use of solid state technology. Infrared radiation with the wavelength of λ=1064 nm is generated by optical pumping of an yttrium aluminum garnet doped with neodymium. The working ultraviolet radiation with the wavelength of $\lambda=213$ nm is formed by cascade transformation with the use of three nonlinear crystals



Efficiency

Thousands of surgeries, including those with the use of personalized ablation technology showed good post-surgery results that were comparable with the published results of the surgeries performed at modern imported excimer laser installations

Stability

The unique system of the laser together with the peculiarities of the 213 nm wavelength ensure obtaining high quality and stable energy from pulse to pulse, from surgery to surgery.

Safety

Absence of hazardous gases in the list of consumables quarantees safety of patients and medical personnel, and may classify the installation as environmentally safe.

Mobility

This is a new opportunity in the world of refractive surgery. The installation may be moved around both within the same building and within mobile health programs. Easy transportation is possible due to small dimensions, low weight and modular design of the installation.

Enhanced surgery possibilities

Customizable software provides to the surgeon a completely new personalized tool for treating glaucoma.

Advantages of the wavelength $\lambda = 213 \text{ nm}$

- The tolerance to hydration of the comea
- Maximum absorption of corneal collagen
- Minimum energy and thermal load on the comea
- Minimum sensitivity to fluctuations in the microclimate of the operating room
- Minimum damaging effect on optical elements

Advantages of Solid State Technology

- No availability of consumable gas mixtures containing toxic fluorine
- No liquid nitrogen to purge the optical path
- · High energy stability for the entire
- High coherence of radiation
- High quality of energy distribution in the spot
- · Environmental, safety and ease of operation



Installation Advantages

- · Easy and easy maintenance
- User-friendly interface
- Multi-level error protection
- Mobility and compactness
- Possibility of transportation by any means of transport
- Long life (at least 10 years)

Clinical Approbation

- Implemented in clinical practice since
- · Successfully works in three clinics of the Russian Federation
- More than 5,000 successful operations were performed

